



Resolution



GP



ENT



PEDIATRIC



TELEMEDICINE

Sync  Vision

iO1 Elite

Digital Otoscope

Safe COVID distance between doctors and patients.





Clinics



Visits



Emergency



Audiology

SyncVision Technology Corp.

3F., No.88, Xingde Rd., Sanchong Dist., New Taipei City 24158, Taiwan
 Tel: +886-2-2278-2222 Fax: +886-2-2278-1111
 E-mail: Info@sync-vision.com Website: www.Sync-Vision.com

Technical features and Specifications

Dimensions	43x40x175mm	File Format	JPEG, MP4
Weight	Approximate 140g	Interface	USB Type C
Display	1.5"full color TFT-LCD	Memory	Micro SD card 32G
Probe Size	Ø3.9mmx 43mm(L)	WIFI standard	802.11
Resolution	Full HD resolution	Power	Built-in rechargeable lithium battery
Focus Type	Fixed focus		3.7V/ 2600mAh / (6 hours runtime)
Light Source	High Luminance White LED		

